



NL 200

## Bluetooth® Neckloop

### FUNCTIONS:

- Compatible with most popular hearing aids (HAC Hearing Aid Compatible to T3 Rating)
- Cordless connection to the mobile phone
- Fully compatible with all Bluetooth® V2.0+ mobile phones
- Range: 10 m
- Integrated lithium-polymer battery
- Talk time: 8 hours
- Standby mode: 2000 hours
- Clear sound, no interference when used with hearing aids using the "T" function
- Integrated microphone
- Volume regulation
- Battery capacity indicator
- Redialling function
- Voice command (depending on mobile phone used)

### CONTENT:

- 1x Neckloop
- 1x Power adapter plug
- 1x User Guide



Warranty : 24 months

EAN No.: 42 60010 593664

Gift box : with EAN Barcode

Export Carton : 6 Gift boxes

Order No : 593664

17.10.2011

Marketing: AUDIOLINE GMBH  
Hellersbergstr. 2a · D-41460 Neuss  
Phone: +49 2131 / 79 90-0  
Sales: +49 2131 / 79 90-10  
Fax: +49 2131 / 27 69 90

